

# Pipelines of Marker Development for Transcriptome-based Exon Capture

## *Part I phylogenomics*

January 17, 2015

Contributors: Sonal Singhal and Ke Bi

For questions or to report bugs, please contact Ke Bi ([kebi@berkeley.edu](mailto:kebi@berkeley.edu))

### Reference:

- [1]. Singhal S. 2013. De novo transcriptomic analyses for non-model organisms: an evaluation of methods across a multi-species data set. *Molecular Ecology Resources* 13:403-416.
- [2]. Bi K, Linderth T, Vanderpool D, Good JM, Nielsen R and Moritz C. 2013. Unlocking the vault: next-generation museum population genomics. *Molecular Ecology* 22:6018-6032.
- [3]. Bi K, Vanderpool D, Singhal S, Linderth T, Moritz C and Good JM. 2012. Transcriptome-based exon capture enables highly cost-effective comparative genomic data collection at moderate evolutionary scales. *BMC Genomics* 13: e403.

The pipelines are deposited in  
<https://github.com/CGRL-QB3-UCBerkeley/MarkerDevelopmentPylogenomics>

---

Scripts included in this pipeline:

[1-PreCleanup](#)

[2-ScrubReads](#)

[3-GenerateAssemblies](#)

[4-AssemblyEvaluation](#)

[5-Annotation](#)

[6-MarkerSelectionTRANS](#)

[6-MarkerSelectionEXONS](#)

\*\*Use "chmod +x script" to make each of these perl scripts executable.

46    \*\*Note: If exon identification is not possible or not desirable, users can use the  
48    entire transcripts for marker development. In this case please use "6-  
MarkerSelectionTRANS". Otherwise please use "6- MarkerSelectionEXONS".

---

50

52 *\*1-PreCleanup\**: Reformats raw cDNA sequencing reads from Illumina HiSeq or  
MiSeq for [2-ScrubReads](#). Specifically, in this step we will remove reads that did not  
54 pass the Illumina quality control filters and modify the sequence identifiers.

56 Dependencies:

FastQC: <http://www.bioinformatics.babraham.ac.uk/projects/fastqc/>

58

### **Input:**

60 Raw sequence data files are grouped and saved in folders named by their sample  
IDs. For instance, three libraries (CGRL\_index1, CGRL\_index15, CGRL\_index40) are  
62 saved under “/home/ke/Desktop/SeqCap/data/rawdata/library/”. Compressed  
fastq sequence files are saved in each of these folders.

64

Fastq files use the following naming scheme:

66 <sample name>\_<barcode sequence>\_L<lane (0-padded to 3 digits)>\_R<read  
number>\_<set number (0-padded to 3 digits)>.fastq.gz

68

For example, in “CGRL\_index15\_CGACCTG\_L006\_R1\_001.fastq.gz”:

70 sample name: CGRL\_index15

barcode sequence: CGACCTG

72 lane (0-padded to 3 digits): 006

read number: 1

74 set number (0-padded to 3 digits): 001

76 #Make a new folder called “raw” under

“~/Desktop/MarkerDevelopment/data/rawdata/”:

78 *ke@NGS:~/Desktop/MarkerDevelopment/data/rawdata\$ mkdir raw*

80 #Copy all these compressed fastq files from each folder (CGRL\_index1,  
CGRL\_index15, CGRL\_index40) to “raw”:

82 *ke@NGS:~/Desktop/MarkerDevelopment/data/rawdata\$ cp*  
*library/CGRL\_index\*/\*.gz raw/*

84

#Check data files in “raw”:

86 *ke@NGS:~/Desktop/MarkerDevelopment/data/rawdata\$ ls raw/\**

*CGRL\_index15\_CGACCTG\_L006\_R1\_001.fastq.gz*

88 *CGRL\_index15\_CGACCTG\_L006\_R2\_001.fastq.gz*

*CGRL\_index1\_TCGCAGG\_L006\_R1\_001.fastq.gz*

90 *CGRL\_index1\_TCGCAGG\_L006\_R2\_001.fastq.gz*

*CGRL\_index40\_TTCGCAA\_L006\_R1\_001.fastq.gz*

92 *CGRL\_index40\_TTCGCAA\_L006\_R2\_001.fastq.gz*

94

### **Commands:**

```

96  #cd to the working directory:
    ke@NGS:~/Desktop/MarkerDevelopment/data/rawdata$ cd ..
98
    #run 1-PreCleanup with fastq evaluation
100 ke@NGS:~/Desktop/MarkerDevelopment/data$ 1-PreCleanup
    ~/Desktop/MarkerDevelopment/data/rawdata/raw/ fastqc
102
    ~/Desktop/MarkerDevelopment/data
104 Output:
    Three new folders will be created under
106 "~/Desktop/MarkerDevelopment/data/rawdata/raw/":
    "pre-clean"
108 "combined"
    "pre-clean/evaluation"
110
    - Folder "pre-clean" contains reformatted raw fastq reads.
112 CGRL_index1_R1.fq
    CGRL_index1_R2.fq
114 CGRL_index15_R1.fq
    CGRL_index15_R2.fq
116 CGRL_index40_R1.fq
    CGRL_index40_R2.fq
118
    - Folder "combined" contains merged, compressed, fastq data files (not used by the
120 following pipeline).
    CGRL_index1_TCGCAGG_L006_R1.fastq.gz
122 CGRL_index1_TCGCAGG_L006_R2.fastq.gz
    CGRL_index15_CGACCTG_L006_R1.fastq.gz
124 CGRL_index15_CGACCTG_L006_R2.fastq.gz
    CGRL_index40_TTCGCAA_L006_R1.fastq.gz
126 CGRL_index40_TTCGCAA_L006_R2.fastq.gz

128 - Folder "evaluation" contains fastQC results for each data file.
    CGRL_index1_R1.fq_fastqc/
130 CGRL_index1_R2.fq_fastqc/
    CGRL_index15_R1.fq_fastqc/
132 CGRL_index15_R2.fq_fastqc/
    CGRL_index40_R1.fq_fastqc/
134 CGRL_index40_R2.fq_fastqc/

```

---

```

136

```

138 \*2-ScrubReads\*: Clean up raw data, which includes trimming for quality, removing  
 140 adapters, merging overlapping reads, removing duplicates and reads sourced from  
 contamination

142 Dependencies:  
 cutadapt: <http://code.google.com/p/cutadapt/>  
 144 COPE: <http://sourceforge.net/projects/coperead/>  
 Bowtie2: <http://sourceforge.net/projects/bowtie-bio/files/bowtie2/>  
 146 FastQC: <http://www.bioinformatics.babraham.ac.uk/projects/fastqc/>  
 FLASh-modified: modified version of FLASh by Filipe G. Vieira.  
 148 <https://github.com/MVZSEQ/Exon-capture>  
 Trimmomatic: <http://www.usadellab.org/cms/?page=trimmomatic>

150 **Input:**

152 1. Reformatted fastq files created by [1-PreCleanup](#):  
 #Check the raw data files:  
 154 *ke@NGS:~/Desktop/MarkerDevelopment/data/rawdata/raw/pre-clean\$ ls \*.fq*  
*CGRL\_index1\_R1.fq*  
 156 *CGRL\_index1\_R2.fq*  
*CGRL\_index15\_R1.fq*  
 158 *CGRL\_index15\_R2.fq*  
*CGRL\_index40\_R1.fq*  
 160 *CGRL\_index40\_R2.fq*

162 2. A fasta file that contains adapter sequences:  
 #Check the format of adapter sequence file:  
 164 *ke@NGS:~/Desktop/SeqCap/denovoTargetCapture/associated\_files \$ less -S*  
*Adapters.fasta*  
 166 *>P7\_index1*  
*CAAGCAGAAGACGGCATACGAGATcctgcgaGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT*  
 168 *>P7\_index2*  
*CAAGCAGAAGACGGCATACGAGATtgcagagGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT*  
 170 *.....*  
*>P5\_index1*  
 172 *AATGATACGGCGACCACCGAGATCTACACcctgcgaACACTCTTCCCTACACGACGCTCTTCCGATCT*  
*>P5\_index2*  
 174 *AATGATACGGCGACCACCGAGATCTACACtgcagagACACTCTTCCCTACACGACGCTCTTCCGATCT*  
 176 *.....*

Note: The header of each adapter sequence has to be named strictly as “**P7\_indexN**”  
 178 or “**P5\_indexN**”. N is the number of index. It is OK to put all adapters in this file but  
 your libraries only use a subset of them.

180 3. Library info file (Tab-delimited txt file):  
 182 #Check the format of Library info file:  
*ke@NGS:~/Desktop/SeqCap/denovoTargetCapture/associated\_files \$ less -S libInfo.txt*  
 184

<i>library</i>	<i>P7</i>	<i>P5</i>
----------------	-----------	-----------

186 *CGRL\_index1* 1  
 187 *CGRL\_index15* 15  
 188 *CGRL\_index40* 40

190 Leave the “P5” column blank if you only have indexes in P7 adapters in the libraries.

192 4. Contaminant file:  
 193 *Escherichia coli* (bacteria + human + other genome resources if desired) genome in  
 194 fasta format.  
 195 This file (e\_coli\_K12.fasta) is saved in  
 196 “~/Desktop/SeqCap/denovoTargetCapture/associated\_files/ecoli/”

198

**Commands:**

200 #Make a new folder called “cleaned\_data” in  
 201 “~/Desktop/MarkerDevelopment/data/”:  
 202 *ke@NGS:~/Desktop/MarkerDevelopment/data\$ mkdir cleaned\_data*

204 #Run [2-ScrubReads](#):  
 205 *ke@NGS:~/Desktop/MarkerDevelopment/data\$ 2-ScrubReads -f*  
 206 *~/Desktop/MarkerDevelopment/data/rawdata/raw/pre-clean/ -o*  
 207 *~/Desktop/MarkerDevelopment/data/cleaned\_data/ -a*  
 208 *~/Desktop/SeqCap/denovoTargetCapture/associated\_files/Adapters.fasta -b*  
 209 *~/Desktop/SeqCap/denovoTargetCapture/associated\_files/libInfo.txt -t*  
 210 */home/ke/Desktop/SeqCap/programs/Trimmomatic-0.32/trimmomatic-0.32.jar -c*  
 211 *~/Desktop/SeqCap/denovoTargetCapture/associated\_files/ecoli/e\_coli\_K12.fasta -e*  
 212 *200 -m 15 -z*

214 Note: I use the default values for most of the arguments. Users should adjust these  
 215 parameters when processing the real datasets.

216

**Output:**

218 1. In “~/Desktop/MarkerDevelopment/data/cleaned\_data/”, six .txt files per  
 219 library are produced:  
 220 For example for library *CGRL\_index1*, the six files are:  
 221 *CGRL\_index1\_1\_final.txt* (left reads)  
 222 *CGRL\_index1\_2\_final.txt* (right reads)  
 223 *CGRL\_index1\_u\_final.txt* (merged or unpaired reads)  
 224 *CGRL\_index1.contam.out* (headers of reads aligned to bacteria)  
 225 *CGRL\_index1.duplicates.out* (headers of duplicated reads)  
 226 *CGRL\_index1.lowComplexity.out* (headers of low complexity reads)

228 2. In “~/Desktop/MarkerDevelopment/data/cleaned\_data/evaluation/”, you can  
 229 find fastQC results for cleaned reads from each library.

230

*\*3-GenerateAssemblies\**: Assemble RNAseq data using Trinity.

Dependencies:

Trinity <http://trinityrnaseq.sourceforge.net>

### **Input:**

For each library, we will concatenate cleaned forward reads (XXX\_1\_final.txt) and unpaired reads (XXX\_u\_final.txt) and name the resulting read data file as XXX\_1\_final.txt.

#Make a new folder called "raw\_assembly" under

"~/Desktop/MarkerDevelopment/data/":

```
ke@NGS:~/Desktop/MarkerDevelopment/data$ mkdir raw_assembly
```

#Concatenate cleaned forward reads and unpaired reads and save them in "raw\_assembly":

```
ke@NGS:~/Desktop/MarkerDevelopment/data$ cat
```

```
cleaned_data/CGRL_index1_1_final.txt cleaned_data/CGRL_index1_u_final.txt | sed 's/\2$/\1/g' > raw_assembly/CGRL_index1_1_final.txt
```

```
ke@NGS:~/Desktop/MarkerDevelopment/data$ cat
```

```
cleaned_data/CGRL_index15_1_final.txt cleaned_data/CGRL_index15_u_final.txt | sed 's/\2$/\1/g' > raw_assembly/CGRL_index15_1_final.txt
```

```
ke@NGS:~/Desktop/MarkerDevelopment/data$ cat
```

```
cleaned_data/CGRL_index40_1_final.txt cleaned_data/CGRL_index40_u_final.txt | sed 's/\2$/\1/g' > raw_assembly/CGRL_index40_1_final.txt
```

#Copy read2 of all libraries to "raw\_assembly"

```
ke@NGS:~/Desktop/MarkerDevelopment/data$ cp
```

```
cleaned_data/CGRL_index*_2_final.txt raw_assembly/
```

### **Commands:**

#Run Trinity on 4 processors.

```
ke@NGS:~/Desktop/MarkerDevelopment/data$ 3-GenerateAssemblies trinity -a
```

```
raw_assembly/ -c 5 -e 4
```

Note: Your laptop may not be able to handle Trinity assemblies.

### **Output:**

There are quite a few intermediate files generated in

"~/Desktop/MarkerDevelopment/data/raw\_assembly/CGRL\_index1/".

"~/Desktop/MarkerDevelopment/data/raw\_assembly/CGRL\_index15/".

"~/Desktop/MarkerDevelopment/data/raw\_assembly/CGRL\_index40/".

#To show final trinity assemblies that are needed for annotation:

```
ke@NGS:~/Desktop/MarkerDevelopment/data$ ls raw_assembly/CGRL_index*/*.fasta
```

```

278 raw_assembly/CGRL_index15/CGRL_index15.fasta
    raw_assembly/CGRL_index1/CGRL_index1.fasta
280 raw_assembly/CGRL_index40/CGRL_index40.fasta

282 #Under "~/Desktop/MarkerDevelopment/data/" make a new folder called
    "annotation" and copy all files shown above to this folder:
284
    ke@NGS:~/Desktop/MarkerDevelopment/data$ mkdir annotation
286 ke@NGS:~/Desktop/MarkerDevelopment/data$ cp
    raw_assembly/CGRL_index*/*.fasta annotation/
288
    #check all files in folder "annotation"
290 ke@NGS:~/Desktop/MarkerDevelopment/data$ ls annotation/*
    annotation/CGRL_index15.fasta
292 annotation/CGRL_index40.fasta
    annotation/CGRL_index1.fasta
294

296 #####

298 When we did step1-3 we used a tiny fraction of the RNAseq data for the
300 purpose of quick demonstration. To better demonstrate how to use the next
    script (4-AssemblyEvaluation) let's sample some more data from each
    individual.
302
    Please do the following before you start working on step 4:
304 ke@NGS:~/Desktop/MarkerDevelopment/data$ cp
    ~/Desktop/MarkerDevelopment/associated_data/CGRL_index*.fasta
306 annotation/

308 #####

310 _____

```



312 *\*4-AssemblyEvaluation\** (Optional): Evaluate the quality of cDNA *de novo* assemblies.  
A few examples of the available functions are shown here.

314

Dependencies:

316 Blat: [http://hgdownload.soe.ucsc.edu/downloads.html#source\\_downloads](http://hgdownload.soe.ucsc.edu/downloads.html#source_downloads)

Blastall:

318 [http://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE\\_TYPE=BlastDocs&DOC\\_TYPE=Download](http://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE_TYPE=BlastDocs&DOC_TYPE=Download)

320

**Input:** Trinity assemblies for all libraries stored in

322 “~/Desktop/MarkerDevelopment/data/annotation/”

324 #Display all items in “~/Desktop/MarkerDevelopment/data/annotation/”

*ke@NGS:~/Desktop/MarkerDevelopment/data/annotation\$ ls*

326

*CGRL\_index15.fasta*

328

*CGRL\_index1.fasta*

*CGRL\_index40.fasta*

330

*a. 4-AssemblyEvaluation BASIC:* function “BASIC” evaluates the quality of in-target  
332 assemblies by reporting basic stats: mean, median, total length, gc%, N50 etc. It also  
generates a distribution of contigs by binned lengths.

334

**Commands:**

336 *ke@NGS:~/Desktop/MarkerDevelopment/data \$ 4-AssemblyEvaluation BASIC -a*  
*annotation/*

338

**Output:**

340 # In folder “~/Desktop/MarkerDevelopment/data/annotation/”, you should get the  
following output files:

342

*CGRL\_index15.hist*

344

*CGRL\_index1.hist*

*CGRL\_index40.hist*

346

*basic\_evaluation.out*

348 **Output:**

1. “XXX.hist” shows distribution of contigs by binned lengths

350

#Display first few lines of the file:

352 *ke@NGS:~/Desktop/MarkerDevelopment/data/annotation\$ head CGRL\_index15.hist*

*200:299 57*

354

*300:399 43*

*400:499 34*

356

*500:599 28*

*600:699 25*

```

358 700:799      23
      800:899      18
360 900:999      24
      1000:1099    13
362 1100:1199    6

364 2. "basic_evaluation.out": results of assembly evaluation
      #Display first few lines of the file:
366 ke@NGS:~/Desktop/MarkerDevelopment/data/annotation$ head
      basic_evaluation.out
368
      b. 4-AssemblyEvaluation ANNOTATABLE: Calculates the percentage of the assembled
370 contigs that get a match in reference. It also calculates average percentage of
      matched bp and mismatches among the matched genes.
372
      Commands:
374 ke@NGS:~/Desktop/MarkerDevelopment/data$ 4-AssemblyEvaluation
      ANNOTATABLE -a annotation/ -b 100 -c
376 ~/Desktop/MarkerDevelopment/associated_data/Xenopus_tropicalis.JGI_4.2.cdna.all.f
      a
378
      Output:
380 #Display results in the output file "annotatable.out":
      ke@NGS:~/Desktop/MarkerDevelopment/data/annotation$ less annotatable.out
382
      Assemblies  total matches(%)  matched bases(%)  avg similarity(%)
384 CGRL_index1    100.00 61.91 78.55
      CGRL_index15 98.00 58.38 77.36
386 CGRL_index40 96.00 68.23 78.17

388 c. 4-AssemblyEvaluation ACCURACY: The percentage of the correctly assembled
      bases estimated using the set of expressed reference transcripts
390
      Commands:
392 ke@NGS:~/Desktop/MarkerDevelopment/data$ 4-AssemblyEvaluation ACCURACY -a
      annotation/ -b 300 -c
394 ~/Desktop/MarkerDevelopment/associated_data/Xenopus_tropicalis.JGI_4.2.pep.all.fa
396
      Output:
      #Display results in the output file "accuracy.out":
398 ke@NGS:~/Desktop/MarkerDevelopment/data/annotation$ less accuracy.out

400 Assemblies  stop codon(%)  gaps(%)
      CGRL_index1 0.000 0.000
402 CGRL_index15 0.692 0.000
      CGRL_index40 0.348 0.000

```

404 d. *4-AssemblyEvaluation CONTIGUITY*: Calculates assembly contiguity (the  
406 percentage of expressed reference transcripts covered by a single, longest  
408 assembled contig) and completeness (the percentage of expressed reference  
transcripts covered by all matched assembled contigs)

410 **Commands:**  
412 *ke@NGS:~/Desktop/MarkerDevelopment/data\$ 4-AssemblyEvaluation CONTIGUITY -*  
414 *a annotation/ -b 300 -c*  
*~/Desktop/MarkerDevelopment/associated\_data/Xenopus\_tropicalis.JGI\_4.2.cdna.all.f*  
*a*

416 **\*\*Note:** that -b in function “CONTIGUITY” refers to the number of randomly selected  
418 sequences from the reference protein database. In functions “BASIC”,  
“ANNOTATABLE” and “ACCURACY” -b refers to the number of randomly selected  
sequences in de novo assemblies\*\*

420 **Output:**  
422 #Display results in the output file “Contiguity.out”:  
424 *ke@NGS:~/Desktop/MarkerDevelopment/data/annotation\$ less Contiguity.out*

<i>Assemblies</i>	<i>complete(%)</i>	<i>contiguity(%)</i>
<i>CGRL_index1</i>	<i>13.46</i>	<i>11.47</i>
<i>CGRL_index15</i>	<i>23.36</i>	<i>23.36</i>
<i>GRL_index40</i>	<i>37.40</i>	<i>30.21</i>

430 \_\_\_\_\_

**\*5-Annotation\*:** annotate assembled contigs using a related reference protein dataset that can be found in Ensembl Genome Browser (<http://www.ensembl.org/index.html>)

Dependencies:  
BLAST+:

[http://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE\\_TYPE=BlastDocs&DOC\\_TYPE=Download](http://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE_TYPE=BlastDocs&DOC_TYPE=Download)

FrameDP: [https://iant.toulouse.inra.fr/FrameDP/cgi-bin/framedp.cgi?\\_wb\\_cfg=/www/iant/FrameDP/cgi-bin/./cfg/FrameDP.cfg&\\_wb\\_session=WBUPAWho&\\_wb\\_main\\_menu=Download&\\_wb\\_function=Download](https://iant.toulouse.inra.fr/FrameDP/cgi-bin/framedp.cgi?_wb_cfg=/www/iant/FrameDP/cgi-bin/./cfg/FrameDP.cfg&_wb_session=WBUPAWho&_wb_main_menu=Download&_wb_function=Download)

exonerate: <http://www.ebi.ac.uk/~guy/exonerate/index.html>

**\*\*Note:** this script works only if you can find a reference database from the EGB. However, if you would like to use NCBI refseq, NR or UniProtKB/Swiss-Prot, modification of this script is needed.


**\*\* Swiss-Prot** (created in 1986) is a high quality manually annotated and non-redundant protein sequence database, which brings together experimental results, computed features and scientific conclusions. UniProtKB/Swiss-Prot is now the reviewed section of the UniProt Knowledgebase.

**\*\*FrameDP:** Sensitive peptide detection on noisy matured sequences. A self-training integrative pipeline for predicting CDS in transcripts which can adapt itself to different levels of sequence qualities.


### Input:

1. download a reference protein dataset from the Ensembl:


Step1. Go to the Ensembl homepage <http://www.ensembl.org/> and click on "Download" located at the top.




Step2. Click on “Download data via FTP” to the left of the download page.




Step3. Select “All” in the “single species data” box in the FTP download page.



Step 4: Find and download the reference. Click on the FASTA link for Protein sequence. In this case we choose *Xenopus tropicalis* as the reference.



Step 5: From FTP server, download reference protein fasta “XXX.pep.all. fa.gz”




## Index of /pub/release-78/fasta/xenopus\_tropicalis/pep/

Name	Size	Date Modified
[parent directory]		
CHECKSUMS	130 B	11/20/14, 4:12:00 PM
README	3.0 kB	11/19/14, 1:16:00 PM
Xenopus_tropicalis.JGI_4.2.pep.abinitio.fa.gz	11.8 MB	11/19/14, 1:16:00 PM
Xenopus_tropicalis.JGI_4.2.pep.all.fa.gz	6.7 MB	11/19/14, 1:02:00 PM

Step 6: unzip the downloaded reference fasta: `gunzip Xenopus_tropicalis.JGI_4.2.pep.all. fa.gz`


Step 7: Find and download the GTF (Gene transfer format (GTF) is a file format used to hold information about gene structure) if there is one available for the reference. In this case we can see that *Xenopus tropicalis* has a GTF so we can download it.




Step 8: unzip the downloaded GTF: `gunzip Xenopus_tropicalis.JGI_4.2.78.gtf.gz`

2. If GTF is not available then you can use Ensembl BioMart tool to obtain a gene annotation file for the reference. For the workshop I will show you how obtain this file from the BioMart tool even though we have downloaded a GTF for the reference.


Step1. Go to the Ensembl homepage <http://www.ensembl.org/> and click on “BioMart” located at the top.



Step2. In the BioMart homepage, select “Ensembl Genes 78” and “Xenopus tropicalis genes (JGI4.2)”.



Step3. Click on “Attributes” icon to the left.



uswest.ensembl.org/biomart/martview/78a82ba2f9f870d4a18592ce7347307e

**Ensembl** BLAST/BLAT | BioMart | Tools | Downloads | Help & Documentation | Blog

New Count Results URL XML Perl

**Dataset**  
Xenopus tropicalis genes (JGI4.2)

**Filters**  
[None selected]  
**Attributes** ←  
Ensembl Gene ID  
Ensembl Transcript ID

**Dataset**  
[None Selected]

Please select columns to be included in the c


☒ **Features** ☐ **Variation**  
☐ **Structures** ☐ **Sequences**  
☐ **Homologs**

☒ **GENE:**

☒ **EXTERNAL:**

☒ **PROTEIN DOMAINS AND FAMILIES:**

542 Step 4. Click on “GENE” to expand the manual. Check on “Ensembl Gene ID” and “Associated Gene Name”.



uswest.ensembl.org/biomart/martview/78a82ba2f9f870d4a18592ce7347307e

**Ensembl** BLAST/BLAT | BioMart | Tools | Downloads | Help & Documentation | Blog

New Count Results URL XML Perl Help

**Dataset**  
Xenopus tropicalis genes (JGI4.2)

**Filters**  
[None selected]  
**Attributes**  
Ensembl Gene ID  
Associated Gene Name

**Dataset**  
[None Selected]

Please select columns to be included in the output and hit 'Results' when ready

☒ **Features** ☐ **Variation**  
☐ **Structures** ☐ **Sequences**  
☐ **Homologs**

☒ **GENE:**

**Ensembl**

☒ Ensembl Gene ID ←

☐ Ensembl Transcript ID

☐ Ensembl Protein ID

☐ Ensembl Exon ID

☐ Description

☐ Chromosome Name

☐ Gene Start (bp)

☐ Gene End (bp)

☐ Strand

☐ Band

☐ Transcript Start (bp)

☐ Transcript End (bp)

☐ Transcription Start Site (TSS)

☐ Transcript length

☐ Transcript Support Level (TSL)

☐ GENCODE basic annotation

☐ APPRIS principal isoform annotation

☒ Associated Gene Name ←

☐ Associated Gene Source

☐ Associated Transcript Name

☐ Associated Transcript Source

☐ Transcript count

☐ % GC content

☐ Gene type

☐ Transcript type

☐ Source (gene)

☐ Source (transcript)

☐ Status (gene)

☐ Status (transcript)


☐ Version (gene)

☐ Version (transcript)


544



546 Step 5. Scroll down the window to find “EXTERNAL”. Click on it to expand the manual. Check on “WikiGene Description”



556 Step6. Click on “Results” icon.



564

566

568

570

Step 7. To export the results, select “CSV” format and check on “Unique results only” box, and then click on “Go”.

uswest.ensembl.org/biomart/martview/78a82ba2f9f870d4a18592ce7347307e

Ensembl BLAST/BLAT | BioMart | Tools | Downloads | Help & Documentation | Blog

Export all results to File  csv ☒ Unique results only

Email notification to

View 10 rows as HTML ☐ Unique results only

Ensembl Gene ID	Associated Gene Name	WikiGene Description
ENSXETG00000002632		
ENSXETG000000026421		
ENSXETG00000008383	golt1b	golgi transport 1B
ENSXETG000000015940		
ENSXETG000000019740		
ENSXETG000000034059		coxsackievirus and adenovirus receptor homolog
ENSXETG000000016254	pafah2	
ENSXETG000000020722	selenbp1	
ENSXETG000000001197	comm7	COMM domain containing 7
ENSXETG000000034332		

Step 8. Save and rename the result to be “*Xenopus.tropicalis*\_gene\_name.txt”. There are three columns, separated by comma:

Ensembl Gene ID, Associated Gene Name, WikiGene Description  
ENSXETG00000008383, golt1b, golgi transport 1B  
ENSXETG000000034059, CARH, coxsackievirus and adenovirus receptor homolog  
ENSXETG00000001197, comm7, COMM domain containing 7

\*\*\*\*\*  
\*\*For this workshop, a reference protein, a GTF and the corresponding biomart gene name file are already downloaded and located in “~/Desktop/MarkerDevelopment/associated\_data/”.

**Input:**

1. A folder that contains all trinity assemblies. These files are located in “~/Desktop/MarkerDevelopment/data/annotation/”
2. Reference protein downloaded from the ensemble:  
Xenopus\_tropicalis.JGI\_4.2.pep.all.fa.
3. Reference biomart gene annotation file:  
Xenopus\_tropicalis\_gene\_name.txt

OR

```

608 4. Reference GTF file:
    Xenopus_tropicalis.JGI_4.2.78.gtf
610
612 Commands:
612 # Run 5-Annotation without a GTF (do not execute the command during the
    workshop, since the runs will take quite a while to finish).
614
614 ke@NGS:~/Desktop/MarkerDevelopment/data$ 5-Annotation -a
616 ~/Desktop/MarkerDevelopment/data/annotation/ -b
616 ~/Desktop/MarkerDevelopment/associated_data/Xenopus_tropicalis.JGI_4.2.pep.all.fa
618 -d ~/Desktop/SeqCap/programs/framedp-1.2.2/ -f
618 ~/Desktop/MarkerDevelopment/associated_data/Xenopus_tropicalis_gene_name.txt -
620 n xenopus -e 1
622
622 ##Copy the annotation results to “~/Desktop/MarkerDevelopment/data”
622 ke@NGS:~/Desktop/MarkerDevelopment/data$ scp -r
624 ~/Desktop/MarkerDevelopment/associated_data/annotation/* annotation/
626
626 Output:
626 For each individual trinity assembly, a new folder is generated under
628 “~/Desktop/MarkerDevelopment/data/annotation/”:
630
630 CGRL_index1_xenopus/
630 CGRL_index14_xenopus /
632 CGRL_index40_xenopus /
634
634 ##The annotated fasta files are named as “XXX_xenopus_annotated.fasta”.
636
636 ke@NGS:~/Desktop/MarkerDevelopment/data/annotation$ ls
636 CGRL_index*/*annotated.fasta
638
638 CGRL_index15_xenopus/CGRL_index15_xenopus_annotated.fasta
640 CGRL_index50_xenopus/CGRL_index50_xenopus_annotated.fasta
640 CGRL_index1_xenopus/CGRL_index1_xenopus_annotated.fasta
642
642 ##make a new folder “probe_design” under
642 “~/Desktop/MarkerDevelopment/data/”.
644 ke@NGS:~/Desktop/MarkerDevelopment/data$ mkdir probe_design
646
646 ##copy all the annotated fasta files to “probe_design”
648 ke@NGS:~/Desktop/MarkerDevelopment/data$ cp
648 annotation/CGRL_index*/*annotated.fasta probe_design/
650
650 ## read and display the first few lines in the annotated fasta file:
652 ke@NGS:~/Desktop/MarkerDevelopment/data/probe_design$ head -4
652 CGRL_index15_xenopus_annotated.fasta

```

```

654 >contig1    gs1_ge432    ENSXETG00000014175    vwa5a NA    5e-57
      TCTCTTACATGGACCCTTCC.....
656 >contig10    5u355_gs356_ge817_3u818 ENSXETG00000004176    mocs2
      molybdenum cofactor synthesis 2 2e-82
658 TGTGCACAGTGTGATGTAG.....

660 For contig1: "gs1" means coding region starts at position 1. "ge432" means coding
      region ends by position 432. No UTRs are present in this contig.
662 "ENSXETG00000014175" is the Ensembl gene ID obtained from Xenopus reference
      database. "vwa5a" is the gene name. "NA" is the wiki gene description and in this
664 case, wiki gene description is missing. "5e-57" is e-value in the BLAST search.

666 For contig10: "5u355" means 5UTR ends by position 355. "gs356" means coding
      region starts at position 356. "ge817" means coding region ends by position 817.
668 "3u818" means 3UTR starts at position 818. "ENSXETG00000004176" is the
      Ensembl gene ID obtained from Xenopus reference database. "mocs2" is the gene
670 name. "molybdenum cofactor synthesis 2" is the wiki gene description. "2e-82" is e-
      value in the BLAST search.
672 ~~~~~
674 Run 5-Annotation with a GTF

676 Commands:
      ke@NGS:~/Desktop/MarkerDevelopment$ 5-Annotation -a
678 ~/Desktop/MarkerDevelopment/data/annotation/ -b
      ~/Desktop/MarkerDevelopment/associated_data/Xenopus_tropicalis.JGI_4.2.pep.all.fa
680 -d ~/Desktop/SeqCap/programs/framedp-1.2.2/ -n xenopus -g
      ~/Desktop/MarkerDevelopment/associated_data/Xenopus_tropicalis.JGI_4.2.78.gtf -e
682 1

684 The output by using GTF is slightly different since the header doesn't have gene
      name descriptions. For example:
686 >contig1    gs1_ge432    ENSXETG00000014175    vwa5a protein_coding    5e-
      57
688 TCTCTTACATGGACCCTTCC.....

690 "gs1" means coding region starts at position 1. "ge432" means coding region ends by
      position 432. No UTRs are present in this contig. "ENSXETG00000014175" is the
692 Ensembl gene ID obtained from Xenopus reference database. "vwa5a" is the gene
      name. "protein_coding" is the type of the gene. "5e-57" is e-value in the BLAST
694 search.

696

698

```

700 **\*6-MarkerSelectionTRANS\***: Find orthologous transcripts in transcriptomes from  
 702 different species and generate input files for probe design. It can be used when exon  
 704 identification is impossible and/or is not preferred.

Dependencies:  
 706 BLAST+:  
[http://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE\\_TYPE=BlastDocs&DOC\\_TYPE=Download](http://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE_TYPE=BlastDocs&DOC_TYPE=Download)  
 708 <http://www.drive5.com/muscle/>  
 710 cd-hit-est: <http://weizhongli-lab.org/cd-hit/>

712 First of all we want to identify orthologous transcripts across transcriptomes from  
 714 different species. We will run the command “6-MarkerSelectionTRANS markers” for  
 this task:

716 **Input:**  
 All annotated transcripts located in  
 718 “~/Desktop/MarkerDevelopment/data/probe\_design”

720 ##  
 ke@NGS:~/Desktop/MarkerDevelopment/data/probe\_design\$ ls  
 722 CGRL\_index15\_xenopus\_annotated.fasta  
 CGRL\_index40\_xenopus\_annotated.fasta  
 724 CGRL\_index1\_xenopus\_annotated.fasta

726 Make a new folder “other\_files” under  
 “~/Desktop/MarkerDevelopment/data/probe\_design/”.  
 728 Use one of the annotated files as a “primary” annotation file. Move the rest to a  
 folder “other\_files”. In the workshop we use CGRL\_index1 as the “primary”  
 730 annotation file.

732 ke@NGS:~/Desktop/MarkerDevelopment/data/probe\_design\$ mkdir other\_files

734 ke@NGS:~/Desktop/MarkerDevelopment/data/probe\_design\$ mv  
 CGRL\_index15\_xenopus\_annotated.fasta CGRL\_index40\_xenopus\_annotated.fasta  
 736 other\_files/

738 **Commands:**  
 # Run 6-MarkerSelectionTRANS markers:  
 740 ke@NGS:~/Desktop/MarkerDevelopment/data\$ 6-MarkerSelectionTRANS markers -f  
 probe\_design/CGRL\_index1\_xenopus\_annotated.fasta -d probe\_design/other\_files/ -a  
 742 1000

744 **Output:**

```

746 #Under "ke@NGS:~/Desktop/MarkerDevelopment/data/probe_design/" a new
    folder called "results" was created by the script.
748 ke@NGS:~/Desktop/MarkerDevelopment/data/probe_design$ cd results/

750 #Markers that passed all filters are stored in "marker_kept.txt". First take a How
    many markers are kept?
752 ke@NGS:~/Desktop/MarkerDevelopment/data/probe_design/results$ wc -l
    marker_kept.txt
754 1050 marker_kept.txt

756 ke@NGS:~/Desktop/MarkerDevelopment/data/probe_design/results$ less -S
    marker_kept.txt
758
Transcript_name: Ensembl Gene ID
760 avgDiv: Average sequence divergence (avg. %mismatches)
    varianceDiv: Variance of sequence divergence
762 avgLength: Average length of the marker
    avgGC: Average CG content of the marker
764 div_CGRL_index15_xenopus_annotated_vs_CGRL_index1_xenopus_annotated:
    sequence divergence between CGRL_index15 and CGRL_index1
766 div_CGRL_index15_xenopus_annotated_vs_CGRL_index40_xenopus_annotated:
    sequence divergence between CGRL_index15 and CGRL_index40
768 div_CGRL_index1_xenopus_annotated_vs_CGRL_index40_xenopus_annotated:
    sequence divergence between CGRL_index1 and CGRL_index40
770

772 #Select the markers that you would like to use for probe design. In this case choose
    the most variable 800 markers and save them in a new file "marker_final.txt"
774
    ke@NGS:~/Desktop/MarkerDevelopment/data/probe_design/results$ tail -800
776 marker_kept.txt > marker_final.txt

778 ++++++
    Now we will use command "6-MarkerSelectionTRANS seq" to generate input fasta
780 files for probe design:

782 Input:
    1. A final set of markers you would like to use for probe design -
784 "~/Desktop/MarkerDevelopment/data/probe_design/results/marker_final.txt"

786 2. A folder containing all trimmed transcripts in fasta format. These files were
    created by "6-MarkerSelectionTRANS markers" and are named as XXX.final2 -
788 "~/Desktop/MarkerDevelopment/data/probe_design/results/"

790
Commands:

```

```

792 # Run 6-MarkerSelectionTRANS seq:
    ke@NGS:~/Desktop/MarkerDevelopment/data$ 6-MarkerSelectionTRANS seq -f
794 probe_design/results/marker_final.txt -d probe_design/results/
    The target size for CGRL_index15_xenopus_annotated.final2 is 700532bp!
796 The target size for CGRL_index1_xenopus_annotated.final2 is 701541bp!
    The target size for CGRL_index40_xenopus_annotated.final2 is 701593bp!
798
    Output:
800 #A new folder "Probe_Design" was created by the script. cd to this folder:
    ke@NGS:~/Desktop/MarkerDevelopment/data/probe_design/results$ cd
802 Probe_Design/

804 #Three fasta sequence files contain sequences of orthologous markers are
    generated and ready for submission for probe design:
806 ke@NGS:~/Desktop/MarkerDevelopment/data/probe_design/results/Probe_Design$
    ls *exonic_targets.txt
808
    CGRL_index15_xenopus_annotated_exonic_targets.txt
810 CGRL_index1_xenopus_annotated_exonic_targets.txt
    CGRL_index40_xenopus_annotated_exonic_targets.txt
812


---


814

```

816 \*6-MarkerSelectionEXONS\*: Find orthologous exons in transcriptomes from different  
species and generate input files for probe design.

818

Dependencies:

820 exonerate: <http://www.ebi.ac.uk/~guy/exonerate/index.html>

cd-hit-est: <http://weizhongli-lab.org/cd-hit/>

822 BLAST+:

[http://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE\\_TYPE=BlastDocs&DOC\\_TYPE=Download](http://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE_TYPE=BlastDocs&DOC_TYPE=Download)

824 MUSCLE: <http://www.drive5.com/muscle/>

826

We will run “6-MarkerSelectionEXONS exons” to identify orthologous exons in  
transcriptomes from each of the species.

830 \*\*If a .gtf file is not available then we will first use a protein and genome reference to  
identify exons from reference species. We will then use the identified exons from the  
832 reference to identify orthologous exons from each of the transcriptomes.

834 \*\*However, if a .gtf file is available then I recommend first run “ParseGTF” to obtain  
exonic sequences from the reference and then run “6-MarkerSelectionEXONS  
836 exons”.

838

+++++

840 First of all we assume no .gtf is available so we have to identify exons using a  
reference protein and reference genome.

842

**Input:**

844 1. Under “~/Desktop/MarkerDevelopment/data/” make a new folder  
“probe\_design\_exons”:

846 *ke@NGS:~/Desktop/MarkerDevelopment/data\$ mkdir probe\_design\_exons/*

848 2. copy all annotated transcripts to  
“~/Desktop/MarkerDevelopment/data/probe\_design\_exons”.

850

*ke@NGS:~/Desktop/MarkerDevelopment/data\$ cp  
852 probe\_design/CGRL\_index1\_xenopus\_annotated.fasta  
probe\_design/other\_files/CGRL\_index\* probe\_design\_exons/*

854

*ke@NGS:~/Desktop/MarkerDevelopment/data\$ cd probe\_design\_exons/*

856

*ke@NGS:~/Desktop/MarkerDevelopment/data/probe\_design\_exons\$ ls*

858 *CGRL\_index15\_xenopus\_annotated.fasta*

*CGRL\_index40\_xenopus\_annotated.fasta*

860 *CGRL\_index1\_xenopus\_annotated.fasta*



862 3. Repeat-masked reference genome  
864 "Xenopus\_tropicalis.JGI\_4.2.dna\_rm.nonchromosomal.fa"

866 4. A reference protein reference "Xenopus\_tropicalis.JGI\_4.2.pep.all.fa";

868 Both 3 and 4 can be downloaded through Ensembl following the instruction above.  
For this workshop these two files are located under  
870 "~/Desktop/MarkerDevelopment/associated\_data".

872 **Command:**  
#Run "6-MarkerSelectionEXONS exons"  
874 ke@NGS:~/Desktop/MarkerDevelopment/data\$ 6-MarkerSelectionEXONS exons -p  
'/home/ke/Desktop/MarkerDevelopment/associated\_data/Xenopus\_tropicalis.JGI\_4.2.  
876 pep.all.fa' -g  
'/home/ke/Desktop/MarkerDevelopment/associated\_data/Xenopus\_tropicalis.JGI\_4.2.  
878 dna\_rm.nonchromosomal.fa' -f  
~/Desktop/MarkerDevelopment/data/probe\_design\_exons -E 1000

880  
\*\* "6-MarkerSelectionEXONS exons" takes very long time to run so please do not run  
882 it during the workshop. Let's skip this step and copy the output files directly from  
"associated\_data":

884  
ke@NGS:~/Desktop/MarkerDevelopment/data\$ cp  
886 ~/Desktop/MarkerDevelopment/associated\_data/probe\_design\_exons/\*.nr  
~/Desktop/MarkerDevelopment/associated\_data/probe\_design\_exons/marker\_\*  
888 probe\_design\_exons/

890

**Output:**  
892 In "~/Desktop/MarkerDevelopment/data/probe\_design\_exons/" there are two  
output files that are relevant for the next step:  
894 1. "marker\_kept.txt": Orthologous exonic markers identified in the three species  
2. "marker\_kept\_one\_exon\_per\_gene.txt" is a subset of "marker\_kept.txt"  
896 ,which contains randomly selected one exon per gene.

898 In both 1 and 2, annotation of each column is explained below:

900 **exon\_name:** Exon ID  
**avgDiv:** Average sequence divergence (avg. %mismatches)  
902 **varianceDiv:** Variance of sequence divergence  
**avgLength:** Average length of the exons  
904 **avgGC:** Average CG content of the exons  
**div\_CGRL\_index15\_xenopus\_annotated\_vs\_CGRL\_index1\_xenopus\_annotated:**  
906 sequence divergence between CGRL\_index15 and CGRL\_index1

**div\_CGRL\_index15\_xenopus\_annotated\_vs\_CGRL\_index40\_xenopus**  
 908 **\_annotated:** sequence divergence between CGRL\_index15 and CGRL\_index40  
**div\_CGRL\_index1\_xenopus\_annotated\_vs\_CGRL\_index40\_xenopus\_annotated:**  
 910 sequence divergence between CGRL\_index1 and CGRL\_index40  
 912  
 914 ++++++  
 Now I will demonstrate how to use *6-MarkerSelectionEXONS* exons when a gtf is  
 916 available.  
 918  
**Command:**  
 920 #Run "ParseGTF":  
 ke@NGS:~/Desktop/MarkerDevelopment/data\$ ParseGTF -f  
 922 ~/Desktop/MarkerDevelopment/associated\_data/Xenopus\_tropicalis.JGI\_4.2.78.gtf -g  
 ~/Desktop/MarkerDevelopment/associated\_data/Xenopus\_tropicalis.JGI\_4.2.dna\_rm.t  
 924 oplevel.fa -o 100 -p 1  
 926 \*\*do not run ParseGTF in the workshop  
 928 **Output:**  
 #Results are stored in "exons.unique" under  
 930 "/home/ke/Desktop/MarkerDevelopment/associated\_data/results"  
 932 #cd to "/home/ke/Desktop/MarkerDevelopment/associated\_data/results"  
 ke@NGS:~/Desktop/MarkerDevelopment/data\$ cd  
 934 ~/Desktop/MarkerDevelopment/associated\_data/results/  
 936 #display at the results  
 ke@NGS:~/Desktop/MarkerDevelopment/associated\_data/results\$ less -S  
 938 exons.unique  
 940 #copy "exons.unique" to  
 "~/Desktop/MarkerDevelopment/data/probe\_design\_exons/"  
 942 ke@NGS:~/Desktop/MarkerDevelopment/associated\_data/results\$ cp exons.unique  
 ~/Desktop/MarkerDevelopment/data/probe\_design\_exons  
 944  
 #copy all annotated transcripts to  
 946 "~/Desktop/MarkerDevelopment/data/probe\_design\_exons/"  
 ke@NGS:~/Desktop/MarkerDevelopment/data\$ cp  
 948 probe\_design/CGRL\_index1\_xenopus\_annotated.fasta  
 probe\_design/other\_files/CGRL\_index\* probe\_design\_exons/  
 950  
 952 **Command:**

```

#run "6-MarkerSelectionEXONS exons" (do not run it in the workshop)
954 ke@NGS:~/Desktop/MarkerDevelopment/data$ 6-MarkerSelectionEXONS exons -p
'/home/ke/Desktop/MarkerDevelopment/associated_data/Xenopus_tropicalis.JGI_4.2.
956 pep.all.fa' -g
'/home/ke/Desktop/MarkerDevelopment/associated_data/Xenopus_tropicalis.JGI_4.2.
958 dna_rm.toplevel.fa' -f ~/Desktop/MarkerDevelopment/data/probe_design_exons -E
1000

```

### **Output:**

962 Same as above:

964 Now we will run command "6-MarkerSelectionEXONS seq" to generate input fasta files for probe design:

966

### **Input:**

968 1. A final set of markers you would like to use for probe design. In this case we choose to use one exon per gene -

970 "~/Desktop/MarkerDevelopment/data/probe\_design/results/  
marker\_kept\_one\_exon\_per\_gene.txt"

972

974 2. A folder containing non-redundant exonic markers in fasta format. These files were created by "6-MarkerSelectionEXONS exons" and are named as XXX\_exon.fa.nr - "~/Desktop/MarkerDevelopment/data/probe\_design\_exons/"

976

### **Commands:**

978 # Run 6-MarkerSelectionEXONS seq:

980 ke@NGS:~/Desktop/MarkerDevelopment/data\$ 6-MarkerSelectionEXONS seq -f  
probe\_design\_exons/marker\_kept\_one\_exon\_per\_gene.txt -d probe\_design\_exons/

982

### **Output:**

984 #A new folder "Probe\_Design" was created by the script. cd to this folder:

986 ke@NGS:~/Desktop/MarkerDevelopment/data/probe\_design\_exons\$ cd  
Probe\_Design/

988 #Three fasta sequence files contain sequences of orthologous exonic markers are generated and ready for submission for probe design:

990

992 ke@NGS:~/Desktop/MarkerDevelopment/data/probe\_design\_exons/Probe\_Design\$ ls  
\*exonic\_targets.txt

994 CGRL\_index15\_xenopus\_annotated\_exonic\_targets.txt

CGRL\_index1\_xenopus\_annotated\_exonic\_targets.txt

996 CGRL\_index40\_xenopus\_annotated\_exonic\_targets.txt